


SimpDo™ Blood Glucose Monitoring System

IVD



- Latest advanced biosensor technology
- Large & clear LCD display
- No wiping/blotting
- Fast & accurate results in 30 seconds
- Low sample volume < 3 µl
- High/low glucose value indicator with alarm
- Easy, 3 step operation

For more details, please visit
www.intecasi.com


SimpDo™ Blood Glucose Monitoring System

Intended Use

The **SimpDo™** blood glucose monitoring system is designed for in vitro diagnostic use only, intended for the quantitative measurement of glucose in fresh capillary whole blood. It is suitable for home or professional use.



1 Insert the test strip to the monitor



2 Apply blood to the test strip



3 The result (in mg/dl or mmol/l) appears in 30 seconds



Specification:

Assay range:	36 ~ 450 mg/dl (2.0~ 25.0 mmol/l)
HCT range:	35% ~ 55%
Relative humidity:	10% ~ 85%
Operating temp range:	15°C ~ 35°C (59°F ~ 95°F)
Memory:	80 results in the memory
Display alarms:	Low below 54 mg/dl (3.0 mmol/l) High over 216 mg/dl (12.0 mmol/l)
Battery:	Approximately 1,000 tests, replaceable battery
Monitor dimensions:	9.8 x 5.0 x 1.5 cm
Weight:	< 50 g (1.76 oz)
Portability:	Test anywhere anytime
Durability:	Auto on/off ensuring longer battery shleflife
Consumables:	Disposable lancets and test strips ensuring safe selftesting

InTec® InTec PRODUCTS, Inc.

Fax: (86)592-6519159 Phone: (86)592-6807173 E-mail: glucose@asintec.com